

Cardio Técnica

Fabrica Argentina de Equipamiento Hospitalario



Designed and manufactured
in Argentina.
Warranty provision of spare parts.

Electrocardiograph

Quality technology for the diagnosis

Electrocardiograph **RG6XX**

Models

Rg603



Recording channels	Printout impresion	Paper width	Display TFT Color	ECG representation in the display	Loading data throught touch keyboard on screen
3	1x12 1x12 + Rhythm 2x6 2x6 +Rhythm 3+2	50mm	4.3"	12 Channels	Yes

Rg606



Recording channels	Printout impresion	Paper width	Display TFT Color	ECG representation in the display	Loading data throught touch keyboard on screen
6	1x12 1x12 + Rhythm 2x6 2x6 + Rhythm 3x4 3x4 + Rhythm 6x2 6 + 3x2	80mm	4.3"	12 Channels	Yes

RG501 Plus



Recording channels	Printout impresion	Paper width	Display LCD	ECG representation in the display	Loading data throught touch keyboard on screen
3	1 1+Ritmo 2 2+Ritmo 3+2	50mm	128X64	3 Channels	No

Principal characteristics

General to all models



High resolution thermal head printer with 8 dot/mm

Digital membrane keyboard

Display TFT color with TouchScreen

Input circuit: Floating

Protect from desfibrillation

Elegant and compact appearance, high mechanical rigidity

Printing of name of professional or institution

Sensitivity: 2.5, 5, 10, 20, 10/5 y auto ganancia

Paper Speed: 5, 10, 12.5, 25 y 50 mm/seg



Operation mode: Manual, Auto and Rhythm

Programmable recording times / selectable from 2 to 10s

Programable filters according to need

AC filter: 50 / 60 Hz

DFT filter : 0.05 / 0.15 / 0.25 / 0.50 Hz

Muscle filter (EMG): 25 / 35 / 45 Hz

Low pass filter: 75 / 100 / 150 Hz

Two - language menu



Key copy of the last ECG recording

Internal memory of studies for later printing

Rs232 port for downloading memory studies to PC



Built-in rechargeable li-ion battery, internal recharger

Large capacity for conducting studies with it's internal battery

AC input 100 to 240 VAC 50/60 Hz



Dimension 270mm x 200mm x 65mm

Weight 1,350 Kgs aproximado without accesories

Weight 1,800kg aproximado whit accesories




Product approved by
Ministry of health and environment
A.N.M.A.T. PM-1230-13.

Diagnostic support



- Perform automatic measurements of times and amplitudes
- Print a presumptive diagnosis
- Printing of the points defined for taking times
- User configurable print formats
- Auto printing of date and time



Accesorios

Includes:

- 1pc Patient Cable
- 1pc Power Cable
- 1pc Roll thermal recording paper
- 4pc (1 set) Limb Electrode
- 6pc (1 set) Chest Electrode
- 1pc Electrocardiograph Bag
- 1pc Quick guide to use



La empresa CardioTécnica SRL tiene certificado su sistema de gestión de calidad por IRAM según norma IRAM-ISO 13485:2016 con número de registro RI-13485-016

Cardio Técnica

Av. Reservas Argentinos 443
C.A.B.A. - Buenos Aires - Argentina
(5411) 4643-9436 / (5411) 4643-9482
+54 9 11 5035 3896

cardiotecnica.com.ar

info@cardiotecnica.com.ar